

Product Datasheet

PPIAL4A antibody (orb1431882)

Description	Rabbit polyclonal antibody to PPIAL4A.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	This PPIAL4A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human PPIAL4A.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	A2BFH1
NCBI	NP_001116540.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.