

Product Datasheet

ACSM4 antibody (orb1432104)

Description	Rabbit polyclonal antibody to ACSM4.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ACSM4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65703
Uniprot ID	POC7M7
NCBI	NP_001073923.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.