

Product Datasheet

SH3BP4 antibody (orb1432440)

Description	Rabbit polyclonal antibody to SH3BP4.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SH3BP4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	107496
Uniprot ID	Q9P0V3
NCBI	NP_055336.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.