

Product Datasheet

JUN (T93) antibody (orb1432907)

Description	Rabbit polyclonal antibody to JUN (T93).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	JUN
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35676
Uniprot ID	P05412
NCBI	NP_002219.1
Dilution Range	IF: 1:10~50, WB: 1:1000
Expiration Date	12 months from date of receipt.