

Product Datasheet

NGFRAP1 antibody (orb1433808)

Description	Rabbit polyclonal antibody to NGFRAP1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	BEX3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	12959
Uniprot ID	Q00994
NCBI	NP_996798.1 , NP_996800.1 , NP_055195.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.