

Product Datasheet

LC3 (APG8A) antibody (orb1433842)

Description	Rabbit polyclonal antibody to LC3 (APG8A).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	MAP1LC3A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14272
Uniprot ID	Q9H492
NCBI	NP_115903.1 , NP_852610.1
Dilution Range	IF: 1:200, WB: 1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.