

Product Datasheet

AADACL2 antibody (orb1434465)

Description	Rabbit polyclonal antibody to AADACL2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	AADACL2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46099
Uniprot ID	Q6P093
NCBI	NP_997248.2
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.