

Product Datasheet

TRIM46 antibody (orb1434691)

Description	Rabbit polyclonal antibody to TRIM46.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TRIM46
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	83424
Uniprot ID	Q7Z4K8
NCBI	NP_001269308.1 , NP_001243530.1 , NP_079334.3 , NP_001243528.1 , NP_001243529.1 , NP_001269307.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.