

Product Datasheet

HRH4 antibody (orb1434839)

Description	Rabbit polyclonal antibody to HRH4.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	HRH4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	44496
Uniprot ID	Q9H3N8
NCBI	NP_001153638.1 , NP_067637.2 , NP_001137300.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.