



Product Datasheet

Nlrp6 antibody (orb1435217)


Description	Rabbit polyclonal antibody to Nlrp6.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	Nlrp6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	97402
Uniprot ID	Q91WS2
NCBI	NP_598707.2
Dilution Range	WB: 1:500, WB: 1:2000, IHC-P-Leica: 1:500, IHC-P-Leica: 1:500
Expiration Date	12 months from date of receipt.



Western blot:All
lanes:Anti-Mouse
Nlrp6 ...



Western blot:Anti-
NLRP6 Antibody at
1:20...



Immunohistochemical
analysis of mouse
sm...