

Product Datasheet

LTB4R2 antibody (orb1435708)

Description	Rabbit polyclonal antibody to LTB4R2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LTB4R2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37942
Uniprot ID	Q9NPC1
NCBI	NP_062813.2 , NP_001158164.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.