



Product Datasheet

ALPL antibody (orb1436162)


Description	Rabbit polyclonal antibody to ALPL.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	ALPL
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57305
Uniprot ID	P05186
NCBI	NP_000469.3 , NP_001120973.2 , NP_001170991.1
Dilution Range	IF: 1:10~50, WB: 1:2000, WB: 1:2000, IHC-P-Leica: 1:500, IHC-P-Leica: 1:500, FC: 1:10~50
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of A...



Western blot:All lanes:Anti-ALPL Antibod...



Western blot:Anti-ALPL Antibody (Center)...