

## Product Datasheet

### Ferritin antibody (orb14366)

**Description**

Ferritin antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested**

ICC, IF, IHC-P, WB

**Applications**
**Immunogen**

KLH conjugated synthetic peptide derived from human Ferritin. Please contact us for the exact immunogen sequence. The peptide is available as orb393956.

**Target**

Ferritin

**Form/Appearance**

10 mM PBS, 0.02% sodium azide

**Concentration**

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

**Storage**

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

Polyclonal antibodies are purified by peptide affinity chromatography

**MW**

20 kDa

**Uniprot ID**

[P29391](#), [P02793](#), [P02792](#)

**NCBI**

[000137](#), [000146](#)

**Entrez**


[2512](#)

**Dilution Range**


IF/ICC: 1:100-800, IHC-P: 1:100-800, WB: 1:200-1000

**Expiration Date**


12 months from date of receipt.



Immunohistochemical staining of paraffin...



IHC-P image of mouse liver tissue using ...



Immunohistochemical staining of paraffin...