

Product Datasheet

TATDN1 antibody (orb1436627)

Description	Rabbit polyclonal antibody to TATDN1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TATDN1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33602
Uniprot ID	Q6P1N9
NCBI	NP_001139632.1 , NP_114415.1
Dilution Range	WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.