



Product Datasheet

CD32b antibody (orb14369)


Description	CD32b antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD32b. Please contact us for the exact immunogen sequence. The peptide is available as orb393959.
Target	CD32b
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	34 kDa
Uniprot ID	Q63203 , P08101
NCBI	001002273
Entrez	2213
Dilution Range	IF/ICC: 1:100-600, IHC-P: 1:100-600, WB: 1:200-1000
Expiration Date	12 months from date of receipt.



IHC-P staining of rat spleen tissue using CD32b antibody.



IF analysis of guinea pig bladder tissue using CD32b antibody.



Western blot analysis of rat urinary bladder tissue using CD32b antibody.