



Product Datasheet

SLC11A2 antibody (orb1437185)


Description	Rabbit polyclonal antibody to SLC11A2.
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	SLC11A2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62266
Uniprot ID	P49281
NCBI	NP_001167596.1, NP_001167599.1, NP_000608.1, NP_001167600.1, NP_001167601.1, NP_001167597.1, NP_001167598.1
Dilution Range	IF: 1:10~50, WB: 1:2000, WB: 1:2000
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of S...



Western blot:All lanes:Anti-SLC11A2 Anti...



Western blot:All lanes:Anti-SLC11A2 Anti...