



Product Datasheet

HSPA5 antibody (orb1437589)


Description	Rabbit polyclonal antibody to HSPA5.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	HSPA5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	72333
Uniprot ID	P11021
NCBI	NP_005338.1
Dilution Range	IF: 1:10~50, WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100, IHC-P: 1:50~100, FC: 1:10~50
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of H...



The anti-HSPA5 Pab is used in Western bl...



Western blot analysis of anti-HSPA5 Pab ...