



Product Datasheet

HtrA1 antibody (orb1437626)


Description	Rabbit polyclonal antibody to HtrA1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	HTRA1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51287
Uniprot ID	Q92743
NCBI	NP_002766.1
Dilution Range	WB: 1:1000, IF: 1:10~50, WB: 1:2000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot:Anti-HtrA1 Antibody (N-term...



Confocal immunofluorescent analysis of H...



Western blot:Anti-HtrA1 Antibody (N-term...