

Product Datasheet

RAP1A antibody (orb1437981)

Description	Rabbit polyclonal antibody to RAP1A.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RAP1A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20987
Uniprot ID	P62834
NCBI	NP_002875.1 , NP_001010935.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.