

Product Datasheet

CT011 antibody (orb1438004)

Description	Rabbit polyclonal antibody to CT011.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GID8
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	26749
Uniprot ID	Q9NWU2
NCBI	NP_060366.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.