

Product Datasheet

ATP6V1B1 antibody (orb1439198)

Description	Rabbit polyclonal antibody to ATP6V1B1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	ATP6V1B1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56833
Uniprot ID	P15313
NCBI	NP_001683.2
Dilution Range	IF: 1:10~50, WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.