

Product Datasheet

TXNDC6 antibody (orb1439653)

Description	Rabbit polyclonal antibody to TXNDC6.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	NME9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36856
Uniprot ID	Q86XW9
NCBI	NP_835231.1
Dilution Range	WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.