

Product Datasheet

FOXL2 antibody (orb1440338)

Description	Rabbit polyclonal antibody to FOXL2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	FOXL2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38772
Uniprot ID	P58012
NCBI	NP_075555.1
Dilution Range	WB: 1:1000, FC: 1:10~50
Expiration Date	12 months from date of receipt.