

Product Datasheet

NTRK1 antibody (orb1440638)

Description	Mouse monoclonal antibody to NTRK1.
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	NTRK1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	87497
Uniprot ID	P04629
NCBI	NP_001012331.1 , NP_001007793.1 , NP_002520.2
Dilution Range	WB: 1:500~1000
Expiration Date	12 months from date of receipt.