

Product Datasheet

ACOX1 antibody (orb1440658)

Description	Mouse monoclonal antibody to ACOX1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	ACOX1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	74424
Uniprot ID	Q15067
NCBI	NP_004026.2 , NP_001171968.1 , NP_009223.2
Dilution Range	IF: 1:10~50, WB: 1:4000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.