

Product Datasheet

SLC9A9 Blocking peptide (orb1441059)

Description	SLC9A9 Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 566-579 of HUMAN SLC9A9
Uniprot ID	Q8IVB4
Expiration Date	12 months from date of receipt.