

Product Datasheet

p27Kip1 (Phospho-T170) Blocking peptide (orb1441146)

Description	p27Kip1 (Phospho-T170) Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 164-176 of HUMAN CDKN1B
Uniprot ID	P46527
Expiration Date	12 months from date of receipt.