

Product Datasheet

PDHE1A (Phospho-S232) Blocking peptide (orb1441250)

Description	PDHE1A (Phospho-S232) Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 226-237 of HUMAN PDHA1
Uniprot ID	P08559
Expiration Date	12 months from date of receipt.