

Product Datasheet

CLEC16A Blocking peptide (orb1442146)

Description	CLEC16A Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 773-785 of HUMAN CLEC16A (HGNC:29013)
Uniprot ID	Q2KHT3
Expiration Date	12 months from date of receipt.