

Product Datasheet

FBXW7 Blocking peptide (orb1442147)

Description	FBXW7 Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 177-188 of HUMAN FBXW7 (HGNC:16712)
Uniprot ID	Q969H0
Expiration Date	12 months from date of receipt.