

Product Datasheet

AFAP1-Y451 Blocking peptide (orb1442693)

Description	AFAP1-Y451 Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 444-457 of HUMAN AFAP1
Uniprot ID	Q8N556
Expiration Date	12 months from date of receipt.