

Product Datasheet

TIMM50 Blocking peptide (orb1442698)

Description	TIMM50 Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 25-36 of HUMAN TIMM50
Uniprot ID	Q3ZCQ8
Expiration Date	12 months from date of receipt.