

Product Datasheet

CASP7(subunit p11) Blocking peptide (orb1442722)

Description	CASP7(subunit p11) Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 208-221 of HUMAN CASP7
Uniprot ID	P55210
Expiration Date	12 months from date of receipt.