

Product Datasheet

MAP1LC3A Blocking peptide (orb1442741)

Description	MAP1LC3A Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 100-121 of HUMAN MAP1LC3A
Uniprot ID	Q9H492
Expiration Date	12 months from date of receipt.