

Product Datasheet

MAN2C1 Blocking peptide (orb1443501)

Description	MAN2C1 Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 128-142 of HUMAN MAN2C1
Uniprot ID	Q9NTJ4
Expiration Date	12 months from date of receipt.