


Product Datasheet


c-Jun (Ab-91) Antibody (orb14436)

Description	c-Jun (Ab-91) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P05412
Expiration Date	12 months from date of receipt.




Immunofluorescence staining of methanol-

...



Immunohistochemical analysis of paraffin...



IHC-P of human breast carcinoma tissue u...