



Product Datasheet

JunB (Ab-79) Antibody (orb14441)

Description	JunB (Ab-79) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P17275
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescent
analysis of H...



Immunohistochemical
analysis of formalin...