



Product Datasheet

JunD (Ab-255) Antibody (orb14443)


Description	JunD (Ab-255) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P17535
Expiration Date	12 months from date of receipt.



Western blot analysis of JunD(Ab-255) an...



Confocal immunofluorescent analysis of H...



Immunohistochemical analysis of formalin...