

Product Datasheet

Guanylyl Cyclase beta 1 (GUCY1B3) Antibody Blocking peptide (orb1445278)

Description	Guanylyl Cyclase beta 1 (GUCY1B3) Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody was selected from the Center region of human GUCY1B3 . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q02153
Expiration Date	12 months from date of receipt.