



Product Datasheet

BCL-2 (Ab-70) Antibody (orb14472)


Description	BCL-2 (Ab-70) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P10415
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Confocal immunofluorescent analysis of H...



Western blot analysis of BCL-2(Ab-70) an...