

Product Datasheet

Bi-APG8b (MAP1LC3B) (Phospho-T93/Y99) Antibody Blocking peptide (orb1449729)

Description	Bi-APG8b (MAP1LC3B) (Phospho-T93/Y99) Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Uniprot ID	Q9GZQ8
Expiration Date	12 months from date of receipt.