

## Product Datasheet

# IGFBP-3(Phospho-S183) Antibody Blocking peptide (orb1449781)

|                          |                                                                                                                                                                                                                                                                         |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | IGFBP-3(Phospho-S183) Antibody<br>Blocking peptide                                                                                                                                                                                                                      |
| <b>Form/Appearance</b>   | Peptides are lyophilized in a solid<br>powder format. Peptides can be<br>reconstituted in solution using the<br>appropriate buffer as needed.                                                                                                                           |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                   |
| <b>Application notes</b> | Synthetic peptide                                                                                                                                                                                                                                                       |
| <b>Protein Sequence</b>  | The synthetic peptide sequence<br>used to generate the antibody was<br>selected from the region of human<br>Phospho-IGFBP-3-pS183. A 10 to<br>100 fold molar excess to antibody is<br>recommended. Precise conditions<br>should be optimized for a particular<br>assay. |
| <b>Uniprot ID</b>        | <b>P17936</b>                                                                                                                                                                                                                                                           |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                                                                                                                                                                                         |