

## Product Datasheet

# KIT(Phospho-Y721) Antibody Blocking peptide (orb1449927)

**Description** KIT(Phospho-Y721) Antibody Blocking peptide

**Form/Appearance** Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

**Note** For research use only

**Application notes** Synthetic peptide

**Protein Sequence** The synthetic peptide sequence used to generate the antibody was selected from the region of human Phospho-KIT-Y721. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.

**Uniprot ID** **P10721**

**Expiration Date** 12 months from date of receipt.