

## Product Datasheet

# MAP2K2(Phospho-T394) Antibody Blocking peptide (orb1449976)

|                          |                                                                                                                                                                                                                                                  |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | MAP2K2(Phospho-T394) Antibody Blocking peptide                                                                                                                                                                                                   |
| <b>Form/Appearance</b>   | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.                                                                                                             |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                            |
| <b>Application notes</b> | Synthetic peptide                                                                                                                                                                                                                                |
| <b>Protein Sequence</b>  | The synthetic peptide sequence used to generate the antibody was selected from the region of human Phospho-MAP2K2-T394. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| <b>Uniprot ID</b>        | <b>P36507</b>                                                                                                                                                                                                                                    |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                                                                                                                                                                  |