

Product Datasheet

AKR1A1 Antibody Blocking peptide (orb1450723)

Description	AKR1A1 Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody was selected from the N-term region of human AKR1A1. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	P14550
Expiration Date	12 months from date of receipt.