

Product Datasheet

HSF1 Sumoylation Site Antibody Blocking peptide (orb1450873)

Description	HSF1 Sumoylation Site Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody is PLVRVKEEPPSPQSPR, containing a predicted sumoylation site from the central region of human HSF1. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q00613
Expiration Date	12 months from date of receipt.