

## Product Datasheet

### LRRIQ4 Antibody Blocking peptide (orb1451926)

|                        |                                                                                                                                      |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | LRRIQ4 Antibody Blocking peptide                                                                                                     |
| <b>Form/Appearance</b> | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| <b>Note</b>            | For research use only                                                                                                                |
| <b>Uniprot ID</b>      | <b>A6NIV6</b>                                                                                                                        |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                      |