

## Product Datasheet

### SDR42E1 Antibody Blocking peptide (orb1454964)

|                        |                                                                                                                                      |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | SDR42E1 Antibody Blocking peptide                                                                                                    |
| <b>Form/Appearance</b> | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| <b>Note</b>            | For research use only                                                                                                                |
| <b>Uniprot ID</b>      | <b>Q8WUS8</b>                                                                                                                        |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                                       |