



Product Datasheet

nNOS (Ab-852) Antibody (orb14557)

Description	nNOS (Ab-852) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P29475
Expiration Date	12 months from date of receipt.



Western blot analysis of nNOS (Ab-852) antibody in rat and mouse brain tissue.



Confocal immunofluorescent analysis of human cells stained with the nNOS (Ab-852) antibody.