

Product Datasheet

CFL2 Antibody Blocking peptide (orb1456408)

Description	CFL2 Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody AP14117a was selected from the N-term region of CFL2. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q9Y281
Expiration Date	12 months from date of receipt.